SUPPLEMENTARY INFORMATION

DMX

Digital Handheld Doppler

For safe use of this device, refer to the Instructions for Use (772637).

772641-02

08/2024


HUNTLEIGH

A MEMBER OF THE ARJO FAMILY

Contents

1. Setup Screens	3
1.1 Vascular and APPG Setup Screens.....	3
1.1.1 Heartrate Display Selection Screen	3
1.1.2 DDNR Hiss Suppression Selection Screen	3
1.2 Vascular Monitoring	4
1.2.1 Vascular Frozen Waveform Options.....	4
1.2.2 Vascular Scroll Options	4
1.2.3 Vascular File Options	4
1.3 APPG Monitoring	5
1.3.1 APPG Live Waveform Screen	5
1.3.2 APPG Trace Only Screen.....	5
1.3.3 APPG Scroll Options	5
1.3.4 APPG File Options	6
1.3.5 APPG Systolic Pressure Adjust Screen	6
1.3.6 APPG Final Review and Save Screen	6
2. Stored file Access	7
2.1 Saved File Confirmation Popup.....	7
2.1.1 Stored File Review Screen.....	7
2.1.2 Demo Files	8
3. Micro SD Card Replacement.....	9

1. Setup Screens


When a Vascular probe is attached, the Setup screen gives access to the Heartrate Display Selection Screen and DDNR Hiss Suppression Screen.

When an APPG probe is attached, the screen gives access to the APPG monitoring functions. ATP and ABI kits are available to purchase as optional extras.

1.1 Vascular and APPG Setup Screens

1.1.1 Heartrate Display Selection Screen

This selection determines if the heartrate figure is displayed when the trace is frozen.


Press  to move the highlight to select display the heartrate  or not display the heartrate .

The sample heartrate figure  appears and disappears accordingly.

Press  to accept your selection or  to discard the changes and return to Setup.

1.1.2 DDNR Hiss Suppression Selection Screen

This selection determines if the DDNR applies hiss suppression to the audio signal.

To disable hiss suppression select .


To enable hiss suppression select .

Press  to move the highlight between selections.

Press  to accept your selection or  to discard the changes and return to Setup.

1.2 Vascular Monitoring

1.2.1 Vascular Frozen Waveform Options


Press  to display the File Options.

Press  to show the scroll options. (This button will only be displayed if more than a screen of data has been captured).

Press  to exit and return to Live Waveform screen.

1.2.2 Vascular Scroll Options

Press  to scroll left / back or  to scroll right / forward through the trace.


Press  to accept the displayed portion of the trace and return to the Frozen Waveform Options screen.

Heart rate is not shown while scrolling.

1.2.3 Vascular File Options

The user can choose to save a waveform or open the folder containing previously stored waveforms for review.

The filename will be generated automatically and a popup temporarily displayed to acknowledge that the file has been saved. Data is stored on the Micro SD card.

Press  to save the **visible portion** of the waveform to permanent storage.


Press  to open the folder containing stored waveforms for review.

Press  to exit and go back to the Frozen Waveform Options screen.

1.3 APPG Monitoring

1.3.1 APPG Live Waveform Screen

The APPG trace runs in real time at the selected time base and is auto scaled to optimise the display height. The time base can be set to 3, 6 or 12 seconds. The currently selected time base is displayed to the right of the chart level with the X axis (3 seconds in the example shown).

Press  to freeze the trace and display the APPG Trace Only screen.

Press  to cycle through the time base options: 3s > 6s > 12s > 3s.

1.3.2 APPG Trace Only Screen

The required time base should be selected before freezing the trace.


Press  to display the File Options.

Press  to show the scroll options. (This button will only be displayed if more than a screen of data has been captured).

Press  to exit the APPG Trace Only screen, return to the APPG Home Screen and restart the trace.

1.3.3 APPG Scroll Options


Press  to scroll left or  to scroll right through the trace.

Press  to accept the displayed portion of the trace and return to the APPG Trace Only screen.

1.3.4 APPG File Options

The user can choose to save a waveform or open the folder containing previously stored waveforms for review.

When saving a waveform, only the portion of the waveform displayed on the screen is actually saved. The user should scroll the desired portion of the waveform into view before saving. The filename will be generated automatically and a popup temporarily displayed to acknowledge that the file has been saved.

Press  to save the visible portion of the waveform to permanent storage

Press  to open the folder containing stored waveforms for review.

Press  to exit back to the APPG Trace Only screen.


1.3.5 APPG Systolic Pressure Adjust Screen

The first pulse detected will be aligned with the vertical line.

If this is not the case, the user can manually adjust the detected systolic pressure using the  or  buttons. This may be useful if the user believes the first detected pulse to be random noise. The trace can be moved in either left or right directions, with the vertical line remaining static. The pressure shown above the vertical line is updated accordingly.

1.3.6 APPG Final Review and Save Screen

The detected or adjusted systolic pressure is displayed to the top left of the screen. The APPG waveform that was saved when the cuff pressure first went above 50 mmHg in the APPG Home screen is displayed.

Press  to save the screen information.

A popup will be displayed momentarily to confirm a successful save .

If data has been saved, press  to exit to the APPG Home Screen.

If the data has not been saved, this following prompt will be displayed.


Press  to exit to APPG Home screen without saving.

Press  to return to Final Review and Save screen.

2. Stored file Access

2.1 Saved File Confirmation Popup

Whenever a file is saved, a popup is temporarily overlaid.

File
Saving

File
Saved


File Save
Failed

Soft keys are unavailable while the popup is on screen.


2.1.1 Stored File Review Screen

A list of stored files is displayed with the most recent file at the top. A bar to the right displays the amount of memory in use.

Press  or  to scroll through the file list.


Press  to select the required file.

Press  to view the selected waveform.

Press  to delete selected waveform.

Press  to confirm delete or  to cancel and return to File options.

Press  or  to scroll through the pages of files.


Press  to return to Live Trace screen.

Highlight  and press  to return to Live Trace screen.

Highlight  and press  to enter the Demo Files screen.

2.1.2 Demo Files

Demo files have been pre-programmed onto the SD Card.

Vascular Demo Files

Press  and  to scroll through the file list and  to select the required file.


Highlight  and press  to return to File Review screen.

Note:

When any Demo file is displayed, the auto switch off timeouts are disabled.

3. Micro SD Card Replacement

1


Disconnect the Doppler from any other equipment.

Use a flat blade screwdriver or similar tool to remove the battery cover.

Remove the batteries.

2

To release the Micro SD card, push it further into the card slot so that it is ejected. Pull the card out of the slot.

Insert the new Micro SD Card into the card slot. Push the card so that it clicks into place.

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Manufactured in the UK by Huntleigh Healthcare Ltd on behalf of;

ArjoHuntleigh AB
Hans Michelsensgatan 10
211 20 Malmö, Sweden

Huntleigh Healthcare Ltd.
35 Portmannor Road, Cardiff, CF24 5HN, United Kingdom
T: +44 (0)29 20485885 sales@huntleigh-diagnostics.co.uk
www.huntleigh-diagnostics.com

Registered No: 942245 England & Wales. Registered Office:
ArjoHuntleigh House, Houghton Hall Business Park, Houghton Regis, Bedfordshire, LU5 5XF
©Huntleigh Healthcare Limited 2021

A Member of the Arjo Family

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.

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Huntleigh Healthcare GmbH