

Lymphoedema
management
at home



HYDROVEN® 12 SERIES

LymphAssist™ Homecare

HUNTLEIGH
A MEMBER OF THE ARJO FAMILY

Manage your lymphoedema treatment more effectively from the comfort of your own home

If you've been diagnosed with primary or secondary Lymphoedema, you will know how important it is to be able to self-manage and continue with your prescribed care plan outside of the traditional hospital and clinical setting.

The LymphAssist™ Homecare compression system from Huntleigh has been designed specifically to compliment your traditional treatment in the comfort of your own home.

The clinically proven and effective LymphAssist™ Homecare system comprises two parts, a simple to use, single key press pump and a range of multi chambered inflatable garments.

The system delivers a soothing and effective treatment pattern which replicates the Manual Lymph Drainage (MLD) technique which would have been previously administered by your therapist.

Clinically proven to be an effective treatment¹



Easy to use

Wear over your normal clothes, and begin treatment with a single button press.



Treatment to fit your lifestyle

With a full cycle lasting only 35 minutes, you can fit your treatment around your busy schedule.



Discreet design

The pump's portable design and QuietConcept™ technology makes home treatment unobtrusive.



Expansion inserts

A range of garment sizes coupled with 'active'² inserts provides an effective treatment for a wide range of limb sizes.

1. Nyree Dunn, Edgar M. Williams, Gina Dolan, and Jane H. Davies. Intermittent Pneumatic Compression for the Treatment of Lower Limb Lymphedema: A Pilot Trial of Sequencing to Mimic Manual Lymphatic Drainage Versus Traditional Graduated Sequential Compression. Lymphatic Research and Biology.
2. Dependent on garment.



Patented treatment profile

The treatment profile replicates the Leduc method of MLD and is made up of two phases; the first phase begins by applying a light pressure up your limb starting at your toes or finger tips, this has the effect of softening the tissue and stimulating the sub surface lymphatic channels.

The second phase then reverses the order by applying a light pressure which clears a pathway closest to the centre of your body. The process then continues down the affected limbs until the pattern has been delivered across the whole length clearing accumulated fluid along the way.

This process takes approximately 35 minutes to deliver.

Arm and upper body lymphoedema treatment designed for you

Providing treatment to the upper body.

Arm

The Lymphassist arm garment uses 12 inflatable chambers to provide treatment from hand to shoulder. Using the principles of manual lymphatic drainage this system gently promotes the healthy movement of excess fluid from the arm.

- Suitable for either arm
- 12 overlapping chambers ensure comfortable treatment and promote movement of lymph in a natural manner
- Easy to put on and take off
- Easy to wipe clean to ensure good hygiene
- 'Active' inserts accommodate larger limbs while contributing to treatment



For bilateral lymphoedema we recommend the use of our jacket garment. However, if preferred, two single arm garments can be attached to a single pump.

Upper body

The Lymphassist Jacket extends treatment beyond a single limb to both arms and the upper body. Excess fluid is moved around the torso to areas where the lymphatic system is operating more effectively.

- Complete bilateral upper body therapy
- Treats both arms and the upper body simultaneously
- 24 overlapping chambers ensure comfortable treatment and promote movement of lymph in a natural manner
- Easy to put on and take off
- Easy to wipe clean to ensure good hygiene
- Front and back chest inserts accommodate different torso shapes
- Arm inserts for larger limbs




Expansion
inserts
available

Upper body garment measuring guide

The following guide is intended as an aid to help the user select the correct size garment. When determining the correct size garment:

- Remove any bandaging, stockings or bulky dressings
- Measure the points as shown on the diagram
- DO NOT pull the measuring tape tight
- Make sure the garment extends beyond the area under treatment. Note that some garments are provided in different lengths
- If circumference measurements exceed those listed, check if inserts are available to provide additional width



*Take the largest measurement when measuring both arms.

LymphAssist (Hydroven 12) Arm Twelve Chamber Arm Garments

Order Code	Type	Arm Length (L)	Upper Arm Circumference (UA)	Hand Circumference (H)
316A68	Full Arm	68cm	30cm to 62cm	15cm to 44cm
316A78	Full Arm	78cm	30cm to 62cm	15cm to 44cm

LymphAssist (Hydroven 12) Arm Insert Twelve Chamber Arm Inserts

Order Code	Type	Arm Length (L)	Additional Width at Upper Arm	Additional Width at Hand
316AI68	Full Arm	68cm	17cm	13cm
316AI78	Full Arm	78cm	17cm	13cm

LymphAssist (Hydroven 12) Jacket 24 Chamber Upper Body Garment

Order Code	Type	Arm Length (L)	Chest Circumference	Upper Arm Circumference (UA)	Hand Circumference (H)
316J78	Long Arm Jacket	68cm	80cm to 90cm	31cm to 60cm	15cm to 41cm
316J68	Short Arm Jacket	58cm	80cm to 90cm	31cm to 60cm	15cm to 41cm

LymphAssist (Hydroven 12) Jacket Inserts (non-inflating)

Order Code	Type	Additional Width	
316JIF	Jacket Insert - Front	10cm	
316JIB	Jacket Insert - Back	10cm	
Order Code	Type	Additional width - Upper Arm	Additional Width - Hand
316JA78	Long Arm Jacket Insert	15cm	5cm
316JA68	Short Arm Jacket Insert	15cm	5cm

Leg and lower body lymphoedema treatment designed for you

Providing treatment to the whole lower body

Leg

The lymphassist leg garment uses 12 inflatable chambers to provide treatment from thigh to ankle. Using the principles of manual lymphatic drainage this system gently promotes the healthy movement of excess fluid from the leg.

- Suitable for either leg
- 12 overlapping chambers ensure comfortable treatment and promote movement of lymph in a natural manner
- Easy to put on and take off
- Easy to wipe clean to ensure good hygiene
- 'Active' inserts accommodate larger limbs while contributing to treatment



For bilateral lymphoedema we recommend the use of our trouser garment. However, if preferred, two single leg garments can be attached to a single pump.

Expansion inserts available

Whole lower body

The lymphassist trousers use the principles of manual lymphatic drainage to apply gentle treatment to both legs, abdomen, hips and genital area.

- Treats both legs and the entire lower body simultaneously
- 24 overlapping chambers ensure comfortable treatment and promote movement of lymph in a natural manner
- Easy to put on and take off
- Easy to wipe clean to ensure good hygiene




Lower body garment measuring guide

The following guide is intended as an aid to help the user select the correct size garment. When determining the correct size garment:

- Remove any bandaging, stockings or bulky dressings
- Measure the points as shown on the diagram above
- DO NOT pull the measuring tape tight
- Make sure the garment extends beyond the area under treatment. Note that some garments are provided in different lengths
- If circumference measurements exceed those listed, check if inserts are available to provide additional width

Lower Body Sizing



*Take the largest measurement when measuring both legs.

LymphAssist (Hydroven 12) Leg

Twelve Chamber Leg Garments

Order Code	Type	Leg Length (L)	Thigh Circumference (T)	Calf Circumference (C)
316L76S	Standard Leg	76cm	55cm to 71cm	33cm to 59cm
316L84S	Standard Leg	84cm	55cm to 71cm	33cm to 59cm
316L76W	Wide Leg	76cm	55cm to 79cm	33cm to 62cm
316L84W	Wide Leg	84cm	55cm to 79cm	33cm to 62cm

LymphAssist (Hydroven 12) Leg Insert

Twelve Chamber Leg Inserts

Order Code	Type	Leg Length (L)	Additional Width at Thigh	Additional Width at Calf
316LI76	Leg Insert	76cm	19cm	14cm
316LI84	Leg Insert	84cm	19cm	14cm

LymphAssist Trousers

24 Chamber Complete Lower Body Garment

Order Code	Type	Leg Length (L)	Waist Circumference	Thigh Circumference	Calf Circumference
316T78	Trousers	98cm (abdomen section additional 29cm)	110cm to 120cm	50cm to 78cm	33cm to 62cm

LymphAssist Trouser Inserts

Order Code	Type	Leg Length (L)	Additional Width at Waist	Additional Width at Thigh	Additional Width at Calf
316TI	Trousers Insert - Single Leg	98cm (abdomen section additional 29cm)	17cm	17cm	17cm
316TIS	Trouser Insert Set (2 x leg inserts)		34cm	17cm	17cm



"LymphAssist has made an incredible difference and improvement to my condition.

After my first session I realised how heavy my limbs had felt before. I didn't expect the liberating feeling it would give me after just one session – it was quite emotional.

I now own a LymphAssist and use it every day. Prior to this my treatment involved compression garments and daily exercise. I still maintain this but using LymphAssist before exercise has made it so much easier."

Mark - LymphAssist™ user

"Having LymphAssist at home is an example of value base health care, enabling patients to take control of their lymphoedema self-management within the comfort of their home at a time and frequency of their choosing.

LymphAssist can thus reduce the need for hospital appointments, as well as travel costs and time off work."

Karen Morgan RGN Bsc (hons) PGCEd NMP
National Lymphoedema Education
and Research Lead in Wales
Lymphoedema Network Wales

As a proud member of the Arjo family, we have been committed to supporting healthcare professionals in improving outcomes and enhancing patient wellbeing since 1979. We do this through our proven solutions for Vascular Assessment & Treatment and Fetal & Patient Monitoring. With innovation and customer satisfaction as our guiding principles, we strive for clinical excellence and improved performance, for life.

Manufactured and distributed by Huntleigh Healthcare Ltd.
35 Portmanmoor Road, Cardiff, CF24 5HN, United Kingdom
T: +44 (0)29 20485885 sales@huntleigh-diagnostics.co.uk
www.huntleigh-diagnostics.com

Registered No: 942245 England & Wales. Registered Office:
ArjoHuntleigh House, Houghton Hall Business Park, Houghton Regis, Bedfordshire, LU5 5XF
©Huntleigh Healthcare Limited 2019

A Member of the Arjo Family

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.

AW: 1001043-4