# HUNTLEIGH

## MD200

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Mode d'emploi

Bruksanvisning

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hwijzing

# INSTRUCTIONS FOR USE

Bruksan

atları

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

## WARNING/CAUTIONS & SAFETY

**Do not** use Dopplers in the presence of flammable gases such as anaesthetic agents.

**Do not** use in the sterile field unless additional barrier precautions are taken.

### Do Not :

- Immerse in any liquid,
- use solvent cleaner,
- use high temperature sterilising processes (such as autoclaving),
- use E-beam or gamma radiation sterilisation.

**Do not** use in the eye.

**Do not** use vascular probes to detect the fetal heart.

**WARNING:** To avoid the risk of electric shock, this equipment must only be connected to supply mains with protective earth.

Avoid positioning the equipment such that access to the IEC power connector becomes difficult.

Any equipment connected to the serial port should comply with EN60601-1, EN60950, EN60065, EN60335 or EN61010.

Dopplex Dopplers are screening tools to aid the healthcare professional and should not be used in place of normal vascular or fetal monitoring. If there is doubt as to vascularity or fetal well-being after using the unit, further investigations should be undertaken immediately using alternative techniques.

We recommend that exposure to ultrasound should be kept As Low As Reasonably Achievable - (ALARA guidelines). This is considered to be good practice and should be observed at all times.

## GENERAL OPERATION & MAINTENANCE

Refer to diagram on inside front cover for Doppler Measuring sites and Recommended Probes.

### During Use

The product will automatically switch off after 5 minutes of no signal, except when connected to Dopplex Reporter software.

When the battery symbol appears on LCD display, the battery needs recharging.

MD200 can be powered by the mains when the battery is low. A standard IEC type socket is mounted on the rear panel.

The internal batteries cannot be overcharged.

### After Use

1. Press and release the On/Off button. If you forget to switch the unit off, it will automatically shut-off.
2. Refer to the cleaning section before storing or using the unit on another patient.
3. Probes and gel can be stored in the **MD200**.

The control unit and the body of the probe are robust and require no special handling. However, the probe tip is delicate and must be handled with care.

### Cleaning

Remove excess gel before parking the probe. Clean with a damp cloth impregnated with mild detergent. Do not allow fluid to seep into either unit.

For disinfection use a soft cloth with Sodium Hypochlorite 1000ppm or alcohol.

Alcohol saturated swabs can be used on the probe tip.

Please be sure to check your local control of infection policies, or any equipment cleaning procedures.

Phenolic, detergent based disinfectants containing cationic surfactants, ammonia based compounds, or antiseptic solutions such as Steriscol or Hibiscrub should never be used on any part of the system.

### Coupling Gel

Use water based ultrasound gel **ONLY**.

## OPERATING INSTRUCTIONS VASCULAR

The **MD200** will select vascular mode when a vascular probe is connected to the control unit. In this mode, bi-directional blood flow rate and direction are indicated on bar graphs (4 levels in each direction) and blood flow is audible in the loudspeaker. Flow towards the probe is in the left speaker and flow away from the probe is in the right speaker. The same applies to the stereo headphones when used. Probe frequency is displayed together with the bar graphs.

### Clinical Use

Apply a liberal amount of gel on the site to be examined. Place the probe at 45° to the skin surface over the vessel to be examined. Adjust the position of the probe to obtain the loudest audio signal. High pitched pulsatile sounds are emitted from arteries while veins emit a non-pulsatile sound similar to a rushing wind.

For best results, keep the probe as still as possible once the optimum position has been found. Adjust the audio volume as required.


### Waveform Recording

Separated waveform outputs are provided in analogue and digital formats. Analogue signals are provided for connection to a single channel chart recorder giving a combined bi-directional waveform, or to a dual channel recorder for separate forward and reverse flow waveforms. (Pin-out details on request).

Digital signals are provided via the RS232 port for printing separate waveforms on a **Dopplex Printa II** or for communicating with **Dopplex Reporter** software package.

### Cal Function

The baseline and sensitivity of the chart recorder can be set up using the Cal function. This generates a zero velocity baseline and a sequence of bi-directional pulses as shown below:

## OPERATING INSTRUCTIONS VASCULAR

### Gain Control

To cope with the wide variety of signals detectable using your **MD200**, a gain control allows you to optimise the bar-graph display and increase the height of the waveforms.

Using the **Mode** button on the **MD200**, gain can be set to x1, x2, x4, x8.

### Connection to Printa II

Hard copy printing is automatically selected when the plug of the interface buffer box is inserted into the RS232 socket on the top panel of the **MD200**. Printing is then initiated by using the **Start/Stop** button.

*Note: Do not switch off the MD200 during printing as this will erase any information stored in the buffer box.*

### Connection to Software

The **MD200** can be connected to the Dopplex Reporter Software Package, via the RS232 interface. The Dopplex Reporter software package is available as an accessory. Communicating is then initiated by using the **Start/Stop** button.

### Probes

Six probes are available for vascular examinations:

<b>VP4HS</b>	4MHz $\pm 1\%$ for deep lying vessels
<b>VP5HS</b>	5MHz $\pm 1\%$ for deep lying vessels and oedematous limbs
<b>VP8HS</b>	8MHz $\pm 1\%$ for peripheral vessels
<b>VP10HS</b>	10MHz $\pm 1\%$ for specialist superficial applications.
<b>EZ8</b>	8MHz $\pm 1\%$ "Widebeam" for peripheral vessels.
<b>IOP8</b>	8MHz $\pm 1\%$ for Intraoperative use.

## OPERATING INSTRUCTIONS OBSTETRICS

Obstetric mode is automatically selected when an obstetric probe (OP2HS/OP3HS) is connected. In this mode, **MD200** provides fetal heart rate (FHR) display with 3 operating modes, and outputs FHR for printing on the Printa II.

### Clinical Use

Apply a liberal amount of gel to the abdomen. Place the faceplate of the probe flat against the abdomen above the symphysis pubis. Adjust the probe to obtain an optimum audio signal ideally by angling the probe around. Avoid sliding it over the skin.

In early pregnancy a full bladder may improve sound detection. In later pregnancy the best signals are generally located higher on the abdomen. The fetal heart sounds like a galloping horse at approximately twice the maternal rate. A wind-like sound is heard from the placenta.

### Standard Mode

In this mode the FHR, averaged over 4-heart beats, is displayed on the 3-digit readout. The LCD displays an outline heart symbol.

### Smoothed Mode

This mode is used to obtain more stable heart rate readings. In this mode, FHR is averaged over 8 beats. The LCD displays a solid heart symbol.

### Manual Mode

This mode is used when a fetal heart beat is audible in the loudspeaker or headphones but, due to noise or a low signal level, the **MD200** cannot reliably calculate the heart rate. In this mode, the heart rate can be manually counted over a period of 10 audible heart beats (see below). The **MD200** will automatically calculate and display the derived FHR on the LCD. The LCD displays a clock symbol.

## OPERATING INSTRUCTIONS OBSTETRICS

### Mode Selection

### Use of Manual Mode

1. Press and hold **Start/Stop** button and immediately count the audible heart beats, counting the first beat as the button is pressed. The LCD displays the flashing clock symbol and the FHR reading is shown as three dashes.
2. Release the Start/Stop button immediately on the count of 10 (i.e. After nine beat intervals). The **MD200** will automatically calculate the derived FHR averaged over the 10 beat period and display the result. This rate value is retained until the measurement is repeated or the unit is switched off. If the button is held for a period less than about 3 seconds the display will clear the previous rate value and reset.

### Probes

OP2HS	2MHz $\pm 1\%$
OP3HS	3MHz $\pm 1\%$

## WARRANTY & SERVICE

Huntleigh Healthcare's standard terms and conditions apply to all sales. A copy is available on request.

These contain full details of warranty terms and do not limit the statutory rights of the consumer.

### Service Returns

If for any reason your Dopplex Unit is being returned, please:

1. Clean the product, as described in the cleaning section.
2. Pack it in suitable packing.
3. Attach the decontamination certificate (or other written statement declaring that the product has been cleaned) to the outside of the package.
4. Mark the package "Service Department - MD200"

### UK

For further details refer to the NHS document HSG(93) 26.

For service, maintenance and any questions regarding this, or any other Huntleigh Healthcare Dopplex product, please contact:

Huntleigh Healthcare Ltd  
Diagnostic Products Division  
35 Portmanmoor Road  
Cardiff  
CF24 5HN UK  
Tel : +44 (0)29 20485885  
Fax: +44 (0)29 20492520

Or your local distributor.

*Manufactured in the UK by Huntleigh Healthcare Ltd. As part of the ongoing development programme the company reserves the right to modify specifications and materials of the MD200 without notice.*

*Dopplex, Huntleigh and 'H' logo are registered trademarks of Huntleigh Technology Ltd 2004.*

©Huntleigh Healthcare Ltd. 2004

## END OF LIFE DISPOSAL

This symbol signifies that this product, including its accessories and consumables is subject to the WEEE (Waste Electrical and Electronic Equipment) regulations and should be disposed of responsibly in accordance with local procedures.

## Electromagnetic Compatibility

Make sure the environment in which the Doppler is installed is not subject to strong sources of electromagnetic interference (e.g. radio transmitters, mobile phones). This equipment generates and uses radio frequency energy. If not installed and used properly, in strict accordance with the manufacturer's instructions, it may cause or be subject to interference. Type-tested in a fully configured system, complies with IEC 60601- 1-2, the standard intended to provide reasonable protection against such interference.

Whether the equipment causes interference may be determined by turning the equipment off and on. If it does cause or is affected by interference, one or more of the following measures may correct the interference:

- Reorienting the equipment
- Relocating the equipment with respect to the source of interference
- Moving the equipment away from the device with which it is interfering
- Plugging the equipment into a different outlet so that the devices are on different branch circuits

## TECHNICAL SPECIFICATION

<b>Product Name:</b>	<b>MD200 Desktop Doppler</b>		
<b>Model No.:</b>	<b>MD200</b>		
<b>Physical Characteristics:</b>	<b>Weight :</b> 1448 gms	<b>Width :</b> 227mm	<b>Depth :</b> 79mm
	<b>Height :</b> 240mm		
<b>Audio Output:</b>	<b>Speaker:</b>	500mW rms typical	
	<b>Headphones :</b>	25 mW rms (32Ω headphones) (max applied voltage +9Vdc)	
<b>RS232:</b>		8 pin sub miniature DIN socket, provides analogue and digital waveform outputs (max Applied voltage +5Vdc).	
<b>Analogue:</b>	Zero crosser, 3.5V full scale per channel (forward and reverse). Conversion factor automatically adjusted to give full scale outputs at $\pm 0.1\%$ ( $\pm 10\%$ ) of probe frequency (e.g. $\pm 8\text{kHz}$ for VP8HS)		
<b>Digital:</b>	RS232C, Data formatted to interface with Printa II or Recorder software running on compatible computer. Conversion factors automatically set as per analogue outputs.		
<b>CAL Function:</b>	Automatic sequence of CAL levels set to $\pm 0.05\%$ ( $\pm 10\%$ ) and $\pm 0.1\%$ ( $\pm 10\%$ ) of probe frequency (e.g. $\pm 4\text{kHz}$ & $\pm 8\text{kHz}$ for VP8HS) with zerobaseline at start and end of sequence.		
<b>Auto Shut-Off:</b>	After 5 minutes of no signal		
<b>Battery Charging:</b>	Recharge time, approximately 2 hours		
<b>Battery Life:</b>	Typically, 400 x 1 minute examinations from full charge (depending on use)		
<b>Fuse:</b>		T250mAH, 250V	
<b>Power Supply:</b>	European version only - 230V ~ 50Hz 10VA USA version only - 120V ~ 60Hz 10VA Saudi Arabia Version Only - 220V ~ 50/60Hz 10VA ~ : Alternating Current Refer to rear panel for details of power supply option.		
<b>Probe Transmitter Frequencies:</b>	OP2HS :	2.0MHz $\pm 1\%$	
	OP3HS :	3.0MHz $\pm 1\%$	
	VP4HS :	4.0MHz $\pm 1\%$	
	VP5HS :	5.0MHz $\pm 1\%$	
	VP8HS :	8.0MHz $\pm 1\%$	
	VP10HS :	10.0MHz $\pm 1\%$	
	EZ8 :	8.0MHz $\pm 1\%$	
	IOP8:	8.0MHz $\pm 1\%$	
<b>Effective Area of Active Transmitter Element:</b>	OP2HS, OP3HS :	233 mm <sup>2</sup> ( $\pm 15\%$ )	
	VP4HS, VP5HS :	22 mm <sup>2</sup> ( $\pm 15\%$ )	
	VP8HS, VP10HS :	8 mm <sup>2</sup> ( $\pm 20\%$ )	
	EZ8 :	16 mm <sup>2</sup> ( $\pm 20\%$ )	
	IOP8:	12.5 mm <sup>2</sup> ( $\pm 20\%$ )	

## TECHNICAL SPECIFICATION

<b>FHR Performance :</b>	<b>Standard Mode</b> Range - 60-210bpm Resolution - 1bpm Averaging - 4 beats Accuracy - $\pm 3$ bpm	<b>Smoothed Mode</b> Range - 60-210bpm Resolution - 1bpm Averaging - 8 beats Accuracy - $\pm 3$ bpm	<b>Manual Mode</b> Range - 60-210bpm Resolution - 1bpm Averaging - 10 beats Accuracy $\pm 3$ bpm (excluding user error)
--------------------------	---	---	--

**Complies With:** IEC 60601-1:1988/A2:1995, EN 60601-1:1990/A13:1996, UL 60601-1 1st Edition, CAN/CSA C22.2 No. 601.1-M90

WITH RESPECT TO ELECTRIC SHOCK, FIRE AND MECHANICAL HAZARDS ONLY, IN ACCORDANCE WITH UL60601-1 and CAN/CSA C22.2 No 601.1.

<b>Classification:</b>	Type of shock protection	Class 1 mains powered equipment
	Degree of shock protection	Type B equipment 
	Protection against water ingress	Not Protected - IPX0
	Degree of safety in presence of flammable gases	Equipment not suitable for use in presence of flammable gases
	Mode of operation	Continuous

<b>Environmental:</b>	<b>Operation</b>	
	Temperature :	+10°C to +30°C
	Relative Humidity :	30% to 75% (non condensing)
	Pressure :	700 to 1060hPa
	<b>Transport and Storage</b>	
	Temperature :	-10°C to +40°C
	Relative Humidity :	20% to 90% (non condensing)
	Pressure :	700 to 1060hPa

A full technical description is provided in the Service Manual 752333.

Medical Devices Directive 93/42/EEC

**Attention, consult accompanying documents / Instructions to Use**

**Attention: consult this manual. Refer to safety section.**

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® and ™ are trademarks of Huntleigh Technology Limited
Huntleigh Healthcare Ltd. - Diagnostic Products Division
HNE Huntleigh Nesbit Evans Healthcare GmbH
ArjoHuntleigh Inc
© Huntleigh Healthcare Limited 2004