INSTRUCTIONS FOR USE

Hydroven 12

Hydroven 12 LymphAssist

Arm, Leg and Insert Garments

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Instructions for Use · Mode d'emploi · Gebrauchsanleitung · Istruzioni per l'uso · Instrucciones de uso ·
Instruções de uso · Gebruiksaanwijzing · Brugervejledning · Bruksanvisning · Instrukcja użycia · Οδηγίες χρήσης ·
Kullanım için talimatla · Bruksanvisning · Käyttöohjeet ·

Arm, Leg and insert Garments · Atelles de jambe et de bras · Arm- und Beinmanschetten · Tutori per braccia e
gambe · Prendas para brazo y pierna · Acessórios para braços (luva) e pernas (meia) · Arm- en beenmanschetten
· Arm- og bentekstiler · Arm- och benmanschetter · Mankiety na rękę i nogę · Ενδύματα Βραχίονα και Κνήμης ·
Kol ve Bacak Manşonları · Arm og Lår plagg · Arm ja Leg Vaatteet ·

316901-12

06/2022

HUNTLEIGH

A MEMBER OF THE ARJO FAMILY

⚠ To be used only under the direction of a physician.

For use with **Huntleigh Hydroven® 12 / Hydroven® 12 LymphAssist** pump only.
Refer to the **Hydroven 12 / Hydroven® 12 LymphAssist** User Manual for complete information on using the pump and garments.
100% Latex Free.


Garments and Garment Inserts


Arm Garment

Arm Garment with Insert

Leg Garment

Leg Garment with Insert

Applying the Garments

1. Make sure the garment is the correct size for the limb, and that a garment insert is used if necessary.
2. Undo the zip on the garment.
3. If a garment insert is fitted to the garment, fully fasten one of the zips between the garment and insert, leaving the other unfastened.
4. Before applying the garment (and insert, if fitted) to the limb, zip up the first 150 mm (6") of the unfastened garment zip. Put the garment (and insert) onto the limb and fully fasten the zip.
5. Make sure the patient is in a comfortable position with the limb supported or elevated as necessary.
6. Connect the garment(s) to the two push-fit, polarised 12-way connectors on the pump. If only one garment is used, fit the rubber cover to the lower pump connector.
7. Connect the pump to a suitable mains power supply. Switch on the pump.
8. Refer to the Hydroven 12 User Manual for further information on changing the therapy settings on the pump control panel.

⚠ WARNINGS

- Do not apply the garment to the limb unless it is partially zipped, as you may damage the garment zip.
- Do not undo the zips or attempt to remove the garments during the therapy session or you may damage the zips. Make sure the therapy session has stopped and the garments have deflated before you remove the garments.
- Do not stand or walk with the garments fitted.
- In the event of a power failure or fault whereby the garment remains inflated, disconnect the tubeset(s) in order to deflate the garment(s) and then remove the garment(s) from the limb(s).

Cleaning the Garment

Wipe down using a neutral detergent or soap powder at 51°C (120°F).

Do not iron

Do not dry clean

Do not tumble dry

Do not machine wash.
Air dry thoroughly.

Disinfecting the Garment

After cleaning, wipe complete garment over using 70% isopropyl alcohol wipe or a chlorine-releasing agent at 1000 p.p.m. available chlorine.

Do not use phenol or phenol-derivative disinfectant.

Rinse with clean water to remove any residue.
Air dry thoroughly.

Cleaning the Tubeset

Use a soft brush. Air dry only.
Do not immerse in water. Do not machine wash.

This section is only applicable to United Kingdom (UK) market when UK marking is applied to the Arjo medical device labelling.

UK Symbol:

UK marking indicating conformity with UK Medical Devices Regulations 2002 (SI 2002 No 618, as amended).
Figures indicate UK Approval Body supervision.

UK Responsible Person & UK Importer:

Arjo (UK) Ltd., ArjoHuntleigh House, Houghton Regis, LU5 5XF

Is the appointed UK Responsible Person as defined in UK Medical Devices Regulations 2002 (SI 2002 No 618, as amended).

For Northern Ireland (NI) CE marking will still apply until further amendment to applicable regulations.

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©Huntleigh Healthcare Limited 2015
Manufactured for **Huntleigh Healthcare Ltd** on behalf of;
ArjoHuntleigh AB
Huntleigh Healthcare Ltd.