DIOP8

Disposable 8 MHz Intraoperative Doppler Probe

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
Instructions for Use · Mode d'emploi · Gebrauchsanleitung · Istruzioni per l'uso · Instrucciones de uso
Bruksanvisning · Gebruiksaanwijzing · Brugervejledning · Bruksanvisning · Käyttöohjeet · Návod k použití

Disposable 8 MHz Intraoperative Doppler Probe · Sonde Doppler peropératoire à usage unique 8 MHz · Intraoperative 8 Mhz Einweg-Doppler-Sonde · Doppler intraoperatorio da 8 MHz monouso · Sonda Doppler intraoperatoria de 8 MHz desechable · Engångs 8 MHz intraoperativ dopplersond · Wegwerp 8 MHz intraoperatieve Doppler-sonde · Disponibel 8 MHz intraoperativ Doppler-sonde · 8 MHz intraoperativ dopplersonder til engangsbruk · Kertakäyttöinen 8 MHz:n intraoperatiivinen doppler-anturi · Jednorázová intraoperační dopplerovská sonda 8 MHz

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Intended Use and Indications

The DIOP8 surgical Doppler is indicated for use by qualified healthcare practitioners in a sterile condition in the operating theatre for the assessment of vascular blood flow by direct application to the vessel wall.

The DIOP8 is suitable for use on all patient populations.

Contraindications

- Do not use the DIOP8 on the eye.
- Do not use in the presence of flammable gases or oxygen rich environments.
- Do not apply the DIOP8 to the patient when using high frequency (HF) surgical equipment. Make sure when using such equipment that the DIOP8 is not in contact with the patient.

Note. No pressure should be applied to the vessel wall.

Warnings/Cautions & Safety

	The sterile Intraoperative Probes are delicate and should be handled with care. Do not drop or strike against hard surfaces. Avoid excessive tension on the probe cable.
	These probes are supplied STERILE. Always ensure that the packaging is intact and undamaged before opening.
	The DIOP8 is a screening tool to aid the healthcare professional and can not provide a diagnosis. If there is doubt as to the status of the vessel after using the unit, further investigations should be undertaken immediately using alternative techniques.
	We recommend that exposure to ultrasound should be kept As Low As Reasonably Achievable - (ALARA guidelines). This is considered to be good practice and should be observed at all times.
	The DIOP8 is for single use only. DO NOT REUSE. DO NOT STERILISE. Failure to follow this advice may result in failure of the product, compromised sterility and potential harm to the patient.
	This equipment is for use only by suitably qualified healthcare practitioners.
	This equipment must not be modified.

Note. The control unit and adaptor are NOT sterile. Refer to the hand-held Doppler IFU for details of cleaning.

Operation

Connect the probe adaptor to the Dopplex cable, using the alignment mark on the barrel.

Locate the probe adaptor into the probe holder at the side of the Dopplex control unit. The Dopplex control unit can then be mounted onto an IV pole using the clamp as shown in Fig. 1.

Using aseptic procedure, remove the DIOP8 probe from its packaging and insert plug into adaptor. Always maintain sterility of the probe.

Note: The Dopplex control unit will automatically switch off at a predetermined time after switching the unit on. See control unit IFU for more information.

The Dopplex control unit can be switched on by a member of staff positioned outside the sterile field.

Adjust the volume control accordingly.

Ensure that the probe connector is fully engaged in the socket on the probe adaptor.

Correct operation of the system should be confirmed by placing the DIOP8 tip wetted with patient's body fluid in light contact with an artery which is known to have blood flowing through it. A clearly audible pulsatile Doppler sound should be heard. The audibility and amplitude of the signal will be optimised when the probe is at an angle of approximately $45^{\circ} \pm 15^{\circ}$ to the vessel.

The probe can now be used to assess blood flow in other vessels.

Aiming the probe proximally along the line of the vessel at an angle of $45^{\circ} \pm 15^{\circ}$ and ensuring that the tip is fully wetted at all times (Fig.2.) will enable the optimal audio signal to be obtained.

(If you are using a DMX, Multi Dopplex II or Super Dopplex II, an indication of the blood velocity and its direction will be shown on the LCD display.)

Operation

The probe can then be moved along the length of the vessel, noting any change in pitch of the Doppler signal or height of the Doppler waveform. This may be indicative of a change in the lumen area.

The same procedure can then be carried out after a graft has been inserted to confirm that adequate blood flow has been restored. By placing the probe on the vessel distal to the anastomosis, confirmation of distal run off is provided.

Product Labelling

	This symbol signifies that this product, including its accessories and consumables is subject to the WEEE (Waste Electrical and Electronic Equipment) regulations and should be disposed of responsibly in accordance with local procedures.				
	This symbol signifies that this product complies with the essential requirements of the Medical Device Directive (93/42/EEC) - Medical Device Regulation (EU/2017/745)				
	Legal Manufacturer in association with the CE mark in Europe AjoHuntleigh AB, Hans Michelsensgatan 10 211 20 Malmö, Sweden		Shelf Life	3 year from date of manufacture.	
Manufactured by: Huntleigh Healthcare Ltd. 35 Portmanmoo Road, Cardiff, CF24 5HN, United Kingdom T: +44 (0)29 20485885 sales@huntleigh-diagnostics.co.uk www.huntleigh-diagnostics.com					
	Sterile Barrier System (non sterile protective outer layer with an internal sterile barrier)	Rx Only	Federal law restricts this device to sale by, or on the order of a licensed healthcare practitioner.		
	Attention, consult accompanying documents / Instructions for Use			Single patient use	
	Applied part is type CF	Warning		Do not resterilise	
	Medical Device	REF		Batch Code	
	Device Identification code		Cardboard packaging can be recycled.		
	Contents can be recycled		Does not contain Latex		
	Temperature Limitations		Humidity Limitations		

Technical Specifications - IEC60601-1 Classification

Type of shock protection	Internally powered equipment (Doplex control unit)	
Degree of shock protection	Type CF equipment	
Protection against ingress of liquids	Probe Adaptor: Ordinary Equipment DIOP8: Suitable for use in contact with body fluids.	
Degree of safety in presence of flammable gases	Equipment not suitable for use in presence of flammable gases or oxygen rich environments	
Mode of Operation	Continuous	

Ultrasound

P_r	< 1 MPa	
I_{ob}	< 20 mW/cm ²	
$I_{spk, 3}$	< 100 mW/cm ²	
Thermal Indices and Mechanical Index are 0.7 or less for all device settings.		

Standards

Complies With	IEC60601-1
EMC	IEC60601-1-2

Environmental

Operating	
Temperature range	+10 °C to +30 °C
Relative Humidity	10% to 90% (non-condensing)
Pressure	860 to 1060 hPa
Storage	
Temperature range	-10 °C to 40 °C
Relative Humidity	90% maximum
Pressure	860 to 1060 hPa


Fig. 1


Fig. 2

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EN	Warranty and Service Huntleigh Healthcare's standard terms and conditions apply to all sales. A copy is available on request. These contain full details of warranty terms and do not limit the statutory rights of the consumer. Service Returns There are NO USER SERVICEABLE PARTS inside the probe adaptor or probe. If for any reason your DIOP8 or probe adaptor is being returned, please: - Contact Service Dept. to obtain instructions and authorisation for the product to be returned. - Failure to do this may result in the product being returned without investigation. FOR UK For further details refer to the NHS document HSG(93)26. For service, maintenance and any questions regarding this, or any other Huntleigh Healthcare Dopplex product, please contact: Huntleigh Healthcare Ltd, Diagnostic Products Division, 35 Portmanmoor Road, Cardiff, CF24 5HN UK, Tel : +44 (0)29 20485885, Fax: +44 (0)29 20492520 Or your local distributor. FOR USA For service, maintenance and any questions regarding this, or any other Huntleigh Healthcare Dopplex product, please contact: Service Department, ArjoHuntleigh Inc, Suite 250, 2349 West Lake Street, Addison, IL60101, USA T: 1-800-323-1245 option 2 W: www.huntleigh-healthcare.us Or your local distributor.
----	--

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
This section is only applicable to United Kingdom (UK) market when UK marking is applied to the Arjo medical device labelling.

UK Symbol:

UK Responsible Person & UK Importer:

Arjo (UK) Ltd., ArjoHuntleigh House, Houghton Regis. LU5 5XF

Is the appointed UK Responsible Person as defined in UK Medical Devices Regulations 2002 (SI 2002 No 618, as amended).

For Northern Ireland (NI) CE marking will still apply until further amendment to applicable regulations.

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Manufactured in the UK by **Huntleigh Healthcare Ltd** on behalf of;

ArjoHuntleigh AB
Hans Michelsensgatan 10
211 20 Malmö, Sweden

2797

Huntleigh Healthcare Ltd.

35 Portmanmoor Road, Cardiff, CF24 5HN, United Kingdom
T: +44 (0)29 20485885 sales@huntleigh-diagnostics.co.uk
www.huntleigh-diagnostics.com

Registered No: 942245 England & Wales. Registered Office:

ArjoHuntleigh House, Houghton Hall Business Park, Houghton Regis, Bedfordshire, LU5 5XF
©Huntleigh Healthcare Limited 2003-2020

A Member of the Arjo Family

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.

HUNTLEIGH

1001048-4

Distributed in the USA by:

Arjo Inc.
2349 West Lake Street, Suite 250
Addison, IL 60101
T: 800-323-1245
www.huntleigh-healthcare.us
1001057-3

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Huntleigh Healthcare GmbH
Ringsraße 54
47533 Kleve, Germany
Tel: +49 (0) 2821-979 2821
Fax: +49 (0) 2821-979 2822
<https://www.huntleigh.de>
1001058-3

www.huntleigh-diagnostics.com/

www.huntleigh-healthcare.us/